English BreathSensor Cable

Intended Use

The BreathSensor Cable is intended for use during sleep disorder studies as a qualitative measure of respiratory airflow for recording onto a data acquisition system. The BreathSensor Cable is designed specifically for use with Nellcor Puritan Bennett, Inc., Model 970 series BreathSensor Airflow Thermistors. The cable allows interfacing of the BreathSensor Airflow Thermistors with the Recording System.

Caution

US Federal law restricts this device to sale by or on the order of a physician. Connect the cable to an electrically isolated input. Patient injury can occur if the cable is not connected properly.

If the cable is damaged, replace it with a new unit. Do not attempt to repair the cable. Unauthorized repair will damage the product and void the warranty.

WEEE-RoHS/Disposal

If applicable, separate collection for electrical and electronic equipment per EC Directive 2002/96/EC. Otherwise, dispose of this device in accordance with local regulations.

Installation Instructions

The cable is a one piece unit. Choose an appropriately sized BreathSensor to fit the patient. Apply the thermistor to the patient as described in the Directions for Use supplied with the thermistor.

- Insert the plastic tab of the thermistor into the cable connector as far as it will go. The tab fits snugly inside the connector.
- Place the security hole on the thermistor over the male portion of the security snap. See Figure 1.

- Close the snap securely. See Figure 2. Plug the connector into the appropriate jack in the amplifier.

NOTE: BreathSensor Model 976 and 978 Airflow Thermistors have a flat, unfolded tab with circuitry printed on both sides. This tab must remain flat when inserted. Do not fold the tab.

Specifications

Operating Temperature	5°C (40°F) - 40°C (104°F)
Storage Temperature	-20°C (-4°F) - 60°C (140°F)
Operating/Storage Humidity	15 - 95% Non-condensing

Cleaning

Wipe the cable with an isopropyl alcohol pad before first use and between each patient use.

Do not immerse the cable in any liquid. Immersion will damage the product and void the warranty. Steam autoclaving and ethylene oxide gas are not recommended. Steam autoclaving or ethylene oxide gas will damage the product and void the warranty.

Troubleshooting

If you have difficulty using this product, please verify the following:

- The BreathSensor and cable connections are properly attached to both the patient and the recorder.
- Ensure the thermistor elements are not contacting patient skin. The majority of thermistor difficulties involve improper thermistor placement.

Warranty

Pro-Tech Services, Inc. warrants this product to be free from defects in materials and workmanship for a period of six months from the date of sale by Pro-Tech Services, Inc. If the product fails to perform in accordance with this warranty, the sole liability of Pro-Tech Services, Inc. is limited to repair or at its option replacement of the product with no charge for parts or labor. Under no circumstances shall Pro-Tech Services, Inc. be liable for any indirect, special or consequential loss or damage, including economic loss, loss of business or overhead, or any punitive damages, arising from the use of, or the inability to use this product. This warranty is rendered void and Pro-Tech Services, Inc. shall not be liable for conditions resulting from damage, marginal performance or malfunctions caused by misuse, abuse, neglect, improper line voltage, power fluctuations, or any adverse environmental conditions, tampering, unauthorized modifications, adjustments or repairs to the product or its accessories, or any other defect not related to materials or workmanship. This warranty is in lieu of all other warranties, express or implied, and is extended only to the original purchaser. Features and specifications are subject to change without notice.

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