English
CPAP Titration Kit
(Circuit Connectors with Port)

Intended Use

The Pro-Tech CPAP Titration kit is for use with CPAP equipment and polysomnographic recorders. The Circuit Connectors are placed inline with the CPAP breathing circuit and allows connection to a pressure transducer to detect respiratory airflow during CPAP titration.

Caution Statements

US federal law restricts this device to sale by or on the order of a physician.

WEEE-RoHS/Disposal

 If applicable, separate collection for electrical and electronic equipment per EC Directive 2002/96/EC. Otherwise, dispose of this device in accordance with local regulations.

Specifications


Operating Temperature	5°C (40°F) - 37°C (98°F)
Storage Temperature	-20°C (-4°F) - 60°C (140°F)
Operating/Storage Humidity	15 - 95% Non-condensing

Installation

General Guidelines:

Assemble the circuit connectors (1) inline with the CPAP machine and patient mask (orientation is not critical and depends on the tubing in use with the CPAP machine). Put tubing (2) from pressure port closest to the CPAP machine to the negative (-) (4) input of the pressure transducer and the port closest to the CPAP mask to the positive (+) (5) input of the pressure transducer. This configuration will result in the inspiratory wave to go up on the recording unit.

- Connect small bore tubing (2) to the circuit connector pressure ports (1).
- Attach the filter (3) to the pressure transducer.

- Circuit connectors
- Connection tubing with Luer Loks
- Single use filter
- Negative input of pressure transducer
- Positive input of pressure transducer

Cleaning

This CPAP Titration Kit was designed to be reusable and can be sterilized by the same procedure used for the large bore CPAP tubing (it may not be autoclaved). The entire assembly must be completely dry before using again to insure patient safety.

The filters are single use only and need to be changed between each patient use.

Warranty

Pro-Tech Services, Inc. warrants this product to be free from defects in materials and workmanship for a period of 30 days from the date of sale by Pro-Tech Services, Inc. If the product fails to perform in accordance with this warranty, the sole liability of Pro-Tech Services, Inc. is limited to repair or at its option replacement of the product with no charge for parts or labor. Under no circumstances shall Pro-Tech Services, Inc. be liable for any indirect, special or consequential loss or damage, including economic loss, loss of business or overhead, or any punitive damages, arising from the use of, or the inability to use this product. This warranty is rendered void and Pro-Tech Services, Inc. shall not be liable for conditions resulting from damage, marginal performance or malfunctions caused by misuse, abuse, neglect, improper line voltage, power fluctuations, or any adverse environmental conditions, tampering, unauthorized modifications, adjustments or repairs to the product or its accessories, or any other defect not related to materials or workmanship. This warranty is in lieu of all other warranties, expressed or implied and is extended only to the original purchaser. Features and specifications are subject to change without notice.

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