## Thermocouple/Thermistor Airflow Sensor User’s Guide

### Intended Use & Features

Pro-Tech thermocouple and thermistor airflow sensors are intended for use during sleep disorder studies to detect respiratory airflow for recording onto a physiological recorder. The thermocouple sensor plugs directly into the AC coupled electrode jack box, no batteries are required. For the thermistor, refer to the polygraph operator’s manual for guidance in determining proper instrument connection and settings. *Caution - sensors are fragile and must be handled gently!*

### Specifications

Operating Temperature 5°C (40°F) - 30°C (86°F)
Storage Temperature -20°C (-4°F) - 60°C (140°F)
Operating/Storage Humidity 15 - 95% Non-condensing

### Sensor Installation

- General Guidelines** - Sensors are to be installed with the sensor tips exposed to the maximum available airflow. The sensor tips should not touch the skin. Skin contact will significantly attenuate the signal amplitude.

- Cannula Style Airflow Sensors** are installed onto the patient in a manner similar to that of an oxygen cannula. The assembly should be carefully secured using both the sliding cinch tubing and surgical tape. The oral sensor may be bent gently into the desired position.

- Disposable Airflow Sensors** are designed for *single use only* and cannot be adequately cleaned for reuse. The sensor is installed onto the patient in the same manner as the cannula style airflow sensor described above. The disposable airflow sensor plugs into the reusable wire set, which is then plugged into the AC coupled electrode jack box.

- Single Element Airflow Sensors** are installed onto the patient as desired. The assembly should be carefully secured using surgical tape.

- Nosebridge Airflow Sensors** are taped to the bridge of the nose with the sensor tips positioned under the nares for maximum available airflow.

- Nellcor Puritan Bennett ADAM® Shell and Breeze® Airflow Sensors** - Please refer to the Nellcor Puritan Bennett’s operating instructions for proper patient installation and cleaning of the ADAM® Shell and Breeze® CPAP masks. The nasal sensors must be gently positioned into the center of the pillow opening. The oral sensors may be bent gently into the desired positions. ADAM® Shell and Breeze® are registered trademarks of Nellcor Puritan Bennett, Inc.

### Thermocouple Instrument Connections

| Sensor Type                     | BLUE Connector            | WHITE Connector       | YELLOW Connector          | BLACK Connector       |
|---------------------------------|---------------------------|-----------------------|---------------------------|-----------------------|
| <b>One Recording Channel</b>    | Airflow Signal G-1 Input  | Reference G-2         | Not Applicable            | Not Applicable        |
| <b>Two Recording Channels</b>   | Nasal Signal G-1 Input    | Oral Signal G-1 Input | Not Applicable            | Common Ref. G-2 Input |
| <b>Three Recording Channels</b> | #1 Nasal Signal G-1 Input | Oral Signal G-1 Input | #2 Nasal Signal G-1 Input | Common Ref. G-2 Input |

### Thermocouple Instrument Settings

- Sensitivity** - 2 - 30 microvolts/mm. Typical signal amplitudes may vary considerably. Response is dependent on such variables as sensor placement, patient airflow, sleep position, ambient room temperature, CPAP pressure, polygraph characteristics, etc.

- Low Frequency Filter/Time Constant** - .16 Hz or lower (1 second or longer).

- High Frequency Filter** - 15 Hz

### Thermistor Instrument Connections and Settings

Refer to the polygraph operator’s manual.

### Cleaning

These sensors are designed to be reusable and can be cleaned with an appropriate cleaning agent that is not corrosive toward metal or plastic. The entire assembly must be completely dry before using again to insure patient safety.

### Caution Statements

**This product is for diagnostic purposes only and is not to be used as an apnea monitor or in a life supporting or life sustaining situation. Federal law restricts this device to sale by or on the order of a physician.**

### Warranty

Pro-Tech Services, Inc. warrants this product to be free from defects in materials and workmanship for a period of fifty uses *or* six months from the date of purchase, *whichever comes first*. The disposable airflow sensor is warranted for a single use only. Should the sensor fail prematurely, the sole liability of Pro-Tech Services, Inc. is limited to repair or at its option replacement of the product with no charge for parts or labor. Under no circumstances shall Pro-Tech Services, Inc. be liable for any loss or damage, direct, consequential, or incidental, including property damage or personal injury arising from the use of, or the inability to use this product. This warranty is rendered void and Pro-Tech Services, Inc. cannot be held liable for conditions resulting from: damage, marginal performance or malfunctions caused by: misuse, abuse, neglect, improper line voltage, power fluctuations, or any adverse environmental conditions, tampering, unauthorized modifications, adjustments or repairs to the product or its accessories. This warranty is in lieu of all other warranties, expressed or implied and is extended only to the original

purchaser. Features and specifications are subject to change without notice.

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Advena Ltd., PO Box 30, Leominster HR6 0ZQ, UK