English Snore Microphone

INTENDED USE

The snore microphone is intended for use during sleep disorder studies to detect snoring sounds for recording onto a physiological recorder.

CAUTION

US Federal law restricts this device to sale by or on the order of a physician.

WARNINGS

This product is for diagnostic purposes only and is not to be used as an apnea monitor or in a life supporting or life sustaining situation.

Sensor must be connected to an electrically isolated input. Patient injury could occur if not connected properly.

WEEE-RoHS/DISPOSAL

♻ If applicable, separate collection for electrical and electronic equipment per EC Directive 2002/96/EC. Otherwise, dispose of this device in ♻ accordance with local regulations.

SPECIFICATIONS

Operating Temperature 5°C (40°F) - 40°C (104°F)
Storage Temperature - 20°C (-4°F) - 60°C (140°F)
Operating/Storage Humidity 15 - 95% Non-condensing

SENSOR INSTALLATION (PATIENT SETUP)

Secure the snore microphone (A) to the patient’s neck as shown in Figure 1 to obtain the maximum signal output (the Pro-Tech logo should be facing away from the patient). The lead wire (B) should be looped and securely taped in place (C) to provide additional strain relief in the event the patient pulls on the wire.


Figure 1

INSTRUMENT CONNECTIONS

Plug the snore microphone’s connectors into the appropriate input of your recording system.

Note: When disconnecting from the recording system grasp the connectors firmly. Do not pull on the wire leads as this can damage the connectors. Damage caused by pulling on the wire leads is not covered under the warranty.

INSTRUMENT SETTINGS

Sensitivity - 20 to 75µV/mm. Adjustment of the sensitivity up or down is typically required.

Low Frequency Filter / Time Constant - 10 Hz / .015 seconds

High Frequency Filter - 70 Hz.

If these settings do not provide a good signal refer to your system software documentation for specific recommendations.

SAMPLE RECORDER TRACING

Figure 2

CLEANING & DISINFECTING

Wipe clean with a hospital approved disinfectant that is not corrosive to plastic or metal. Snore microphones are not designed to be sterilized and are not liquid tight. Do not immerse the sensor in cleaning solution or allow liquid to get in to the opening (some models). **Immersing the sensor will damage the device and void the product warranty.**

WARRANTY

Pro-Tech Services, Inc. warrants this product to be free from defects in materials and workmanship for a period of 90 days from the date of sale by Pro-Tech Services, Inc. If the product fails to perform in accordance with this warranty, the sole liability of Pro-Tech Services, Inc. is limited to repair or at its option replacement of the product with no charge for parts or labor. Under no circumstances shall Pro-Tech Services, Inc. be liable for any indirect, special or consequential loss or damage, including economic loss, loss of business or overhead, or any punitive damages, arising from the use of, or the inability to use this product. This warranty is rendered void and Pro-Tech Services, Inc. shall not be liable for conditions resulting from damage, marginal performance or malfunctions caused by misuse, abuse, neglect, improper line voltage, power fluctuations, or any adverse environmental conditions, tampering, unauthorized modifications, adjustments or repairs to the product or its accessories, or any other defect not related to materials or workmanship. This warranty is in lieu of all other warranties, expressed or implied and is extended only to the original purchaser. Features and specifications are subject to change without notice.

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All rights reserved.
CE REP Advena Ltd. Hereford HR4 9DQ UK.