**zRIP™ Module  
User Guide**

## **zRIP™ Module User Guide**

Pro-Tech, a Respironics company  
Tel: 800-345-6443, US & Canada  
724-387-4000, International  
Fax: 800-886-0245, US & Canada  
724-387-5012, International

Email:

clinical@respironics.com, US  
Sleepdiagsupintl@respironics.com, International

Manufactured by:  
Pro-Tech Services, Inc  
4338 Harbour Pointe Blvd SW  
Mukilteo, WA 98275 USA

**EC REP** Advena Ltd. Hereford HR4 9DQ UK

## **Intended Use**

The zRIP module is intended for use with Pro-Tech zRIP Respiratory Effort Sensors during sleep disorder studies for the purpose of conditioning and amplifying the thoracic and abdominal effort signals. The Sum channel (*SUM models only*) is the sum of the thoracic and abdominal effort signals and provides a visual aid to help the technologist identify desynchronized patient abdominal and thoracic breathing.

## **Warnings**

The zRIP System is not to be used in the presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide.

To minimize toxic hazards to the environment, dispose of the batteries according to the manufacturer's instructions and in accordance with local regulations. Pro-Tech Services, Inc. is not responsible for the disposal of the batteries.

Do not operate the zRIP System with the battery cover open.

## **Cautions**

This product is for diagnostic purposes only and is not to be used as an apnea monitor or in a life supporting or life sustaining situation. Federal law restricts this device to sale by or on the order of a Physician.

## **Specifications**

Operating Temperature	5°C (40°F) - 40°C (104°F)
Storage Temperature	-20°C (-4°F) - 60°C (140°F)
Operating/Storage Humidity	15 - 95% Non-condensing
Battery Type	AAA (1.5V)
Battery Life	600 hours continuous operation

## **System Components**

Typical zRIP System Components (sold separately or in kits):

- zRIP Effort Sensor, Adult, Large, or
- zRIP Effort Sensor, Adult, or
- zRIP Effort Sensor, Pediatric,
- zRIP Wireset,
- zRIP Module

Two zRIP effort sensors are typically used with each patient, one for the thorax and one for the abdomen. zRIP wiresets are used to connect the zRIP effort sensors to the zRIP module which is then connected to the physiological recorder. However, some physiological recorder systems may accept the zRIP effort sensors without the use of the zRIP module.

## **Battery Check and Replacement**

The zRIP Interface Box uses two AAA batteries. The ready light will show green when the test button is pressed if the batteries are good. If no light comes on when the test button is pressed, replace the batteries.

Note: If the zRIP module is going to be stored for an extend length of time it is recommended that the batteries be removed.

## **System Setup**

### **Sensor Installation**

Adjust effort sensors so they are approximately the size of the patient (thorax and abdomen):

#### **To lengthen the sensor:**

Hold the male buckle, pull the inner layer of the loop down and away from the buckle, then grasp the sensor past the slide lock and pull them apart until the loop is flat;


#### **To shorten the sensor:**

Hold the slide lock and pull the inner layer of the loop toward the male buckle, then grasping the slide lock and the male buckle, pull them apart until the loop is flat.

Snap buckles together in the front. Adjust the sensor so that it is secure, but not overly tight:

To loosen the sensor, grasp the slide lock and lift slightly away from the patient while pulling

towards the loop; To tighten the sensor, grasp the slide lock and pull away from the loop while repeatedly grasping and pulling the outer layer of the loop toward the slide lock.


**Figure 1**

Insure that the sensors are placed appropriately over the abdomen and thorax as shown in Figure 2. Proper placement is essential for the correct presentation of the signals.

**Wireset Connection**

Plug the dual 1.5mm ends of the wireset into the effort sensor connection; plug the single 249 connector into the designated inputs on the zRIP module.

NOTE: The module will power up when abdominal wireset is plugged into the module.

**Figure 2**

**zRIP Module**

The output from the thoracic and abdominal zRIP effort sensors are plugged into the designated inputs of the zRIP module. The zRIP module outputs are plugged into the designated AC coupled inputs of the recording system, i.e., Thoracic, Abdominal, and optional Sum Channel. The blue plug is signal and the white is reference/common.

zRIP Sum: The Sum Channel automatically compensates for effort changes by averaging previous breaths.

zRIP Sum RT: The Sum Channel is the direct sum of the effort signals.

Note: The zRIP modules will not function with wiresets other than those provided specifically for this purpose.

**Instrument Settings**

**Sensitivity** - Approximately 50µv/mm per channel. Adjustment of the sensitivity up or down is typically required. Response is dependent upon such variables as sensor application and patient effort.

**Low Frequency Filter / Time Constant** - .16 Hz / 1 second or longer. Shorter time constants or higher low frequency filter settings will significantly attenuate waveforms. Longer time constants or lower low frequency filter settings will display more information regarding patient effort.

**High Frequency Filter** - 35 Hz

**Cleaning**

Surface cleaning only with water dampened cloth. The module is not designed to be submersed in liquid.

**Warranty**

Pro-Tech Service, Inc. warrants this product to be free of defects in materials and workmanship for a period of one year from the date of purchase. Should any product prove defective in workmanship, materials or performance during the warranty, the sole liability of Pro-Tech Services, Inc. is limited to repair or, at its option, replacement of the product with no charge for parts or labor. Under no circumstances shall Pro-Tech Service, Inc. be liable for any loss or damage, direct, consequential or incidental including property damage or personal injury arising from the use of, or the inability to use this product. This warranty is rendered void and Pro-Tech Service, Inc. cannot be held liable for conditions resulting from: damage, marginal performance or malfunctions caused by: misuse, abuse, neglect, improper line voltage, power fluctuations, or any adverse environmental conditions, tampering, unauthorized modifications, adjustments or repairs to the product or its accessories. This warranty is in lieu of all other warranties, expressed or implied, and is extended only to the original purchaser. Features and specifications are subject to change without notice.

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